

Leaded X-Ray Glass Frame Listing

Rx Safety Optimum frames for leaded lenses and side shields		
Rx Safety Frame		Colour Options
ASW Point Break		Matte Black & Black/Purple
ASW Trifecta		Matte Black & Gunmetal/Blue
CO 360		Crystal Grey
OG 210		Black & Brown
OG 220		Pink/Grey, Purple/White, Black/Clear, Brown/Clear & Clear
OG 225		Navy/Orange, Black/Yellow, Gunmetal/Navy & Gunmetal/Green
UF Cannon		Matte Black
UF Lynx		Black/Teal
UF Sparky		Black/Crystal
UF Twister		Black, Brown Woodgrain & Burgundy
UF Warhead		Matte Black, Brown Tortoiseshell & Grey/Red

Rx Safety frames <u>not</u> recommended (and why)		
ASW Glacier		<i>(7 base)</i>
CO 108		<i>(Shield sits over temple - can't cover)</i>
CO 200		<i>(Shield sits over temple - can't cover)</i>
CO 202		<i>(Shield sits over temple - can't cover)</i>
CO 207		<i>(Shield sits over temple - can't cover)</i>
CO 333		<i>(8 base)</i>
CO 404		<i>(too rigid)</i>
CO 525		
CO IArmour		<i>(8 base)</i>
UF Grenade		<i>(Bad frame, gaps, wide bridge, large step, won't fit into insert)</i>
UF Whirlwind		<i>(Just not quite covers eyeball)</i>